

# TRBP antibody

Catalog Number: 116352

## Product name

TRBP antibody

Recommended Dilution:

WB: 1:500-1:5000

## Specificity

Human, Mouse, Rat; other species not tested.

IP: 1:500-1:5000


IHC: 1:20-1:200

## Antibody description

TRBP Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in HepG2 cells, mouse liver tissue, mouse spleen tissue. Positive IF detected in Hela cells. Positive IHC detected in human brain tissue, human kidney tissue, human ovary tissue, human placenta tissue, human spleen tissue, human testis tissue. Observed molecular weight by Western-blot: 43-46 kDa

IF: 1:10-1:100

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116352(TRBP antibody) at dilution of 1:1000

## Preparation

This antibody was obtained by immunization of TRBP recombinant protein (Accession Number: NM\_134323). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal


## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IF, IP


## Dilutions




Immunohistochemical of paraffin-embedded human brain using Catalog No:116352(TRBP antibody) at dilution of 1:100 (under 40x lens)

# TRBP antibody

Catalog Number: 116352



Immunofluorescent analysis of HeLa cells, using TARBP2 antibody Catalog No:116352 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-TRBP (IP:Catalog No:116352, 4ug; Detection:Catalog No:116352 1:1500) with mouse liver tissue lysate 6400ug.