

Anti-phospho-Nucleophosmin (Thr199) Rabbit antibody



Catalog Number: 167752

Product name

Rabbit IgG

Anti-phospho-Nucleophosmin (Thr199) Rabbit antibody

Applications

WB, IHC-P, ICC/IF, FC, IP

Specificity

Human

Dilutions

WB: 1/2000

Antibody description

phospho-Nucleophosmin (Thr199) Rabbit monoclonal antibody

IHC: 1/50

ICC/IF: 1/50

FC: 1/20

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Thr199 of human Nucleophosmin

IP: 1/20

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of Nucleophosmin (phospho-Thr199) in C6,Ramos using Nucleophosmin (phospho-Thr199) Rabbit mAb(1:1000 diluted)