

ACTRT2 antibody



Catalog Number: 107752

Product name

ACTRT2 antibody

Specificity

Human,Mouse,Rat; other species not tested.

Antibody description

ACTRT2 Rabbit Polyclonal antibody. Positive WB detected in rat testis tissue, mouse testis tissue. Observed molecular weight by Western-blot: 40-42 kDa

Preparation

This antibody was obtained by immunization of ACTRT2 recombinant protein (Accession Number: NM_080431). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



rat testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107752(ACTRT2 Antibody) at dilution of 1:600