

DCLRE1B antibody



Catalog Number: 109752

Product name

DCLRE1B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DCLRE1B Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HepG2 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of DCLRE1B recombinant protein (Accession Number: NM_022836). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:50-1:500

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109752(DCLRE1B antibody) at dilution of 1:800



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog No:109752(DCLRE1B Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)