

GPD1L antibody



Catalog Number: 111062

Product name

GPD1L antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

GPD1L Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, HEK-293 cells, human brain tissue. Observed molecular weight by Western-blot: 38-40kd

Preparation

This antibody was obtained by immunization of GPD1L recombinant protein (Accession Number: NM_015141). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:111062(GPD1L antibody) at dilution of 1:500