

ALOX15 Polyclonal Antibody



Catalog Number: 165062

Product name

ALOX15 Polyclonal Antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

Polyclonal antibody to ALOX15

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human ALOX15 (NP_001131.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - ALOX15 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ALOX15 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - ALOX15 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat lung using ALOX15 antibody at dilution of 1:200



(40x lens).

ALOX15 Polyclonal Antibody



Catalog Number: 165062

Immunohistochemistry - ALOX15 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach cancer using ALOX15 antibody at dilution of 1:200 (40x lens).