

TIP30 antibody



Catalog Number: 116062

Product name

TIP30 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TIP30 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, SW 1990 cells. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of TIP30 recombinant protein (Accession Number: NM_001098520). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116062(TIP30 Antibody) at dilution of 1:600