

TINAGL1 Polyclonal Antibody



Catalog Number: 161462

Product name

TINAGL1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to TINAGL1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 258-467 of human TINAGL1 (NP_071447.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

Validations



Western blot - TINAGL1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TINAGL1 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.