

# OATP14 antibody

Catalog Number: 113462

## Product name

OATP14 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

OATP14 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in PC-3 cells, mouse brain tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 62-70 kDa

## Preparation

This antibody was obtained by immunization of OATP14 recombinant protein (Accession Number: NM\_017435). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:113462(OATP14 antibody) at dilution of 1:1000



IP Result of anti-OATP14 (IP:Catalog No:113462, 4ug; Detection:Catalog No:113462 1:800) with mouse brain tissue lysate 4000ug.