

Phospho-Syk-Y352/ZAP70-Y319 pAb



Catalog Number: 166462

Product name

Phospho-Syk-Y352/ZAP70-Y319 pAb

Specificity

Human

Antibody description

Polyclonal antibody to Phospho-Syk-Y352/ZAP70-Y319

Preparation

Antigen: A synthetic phosphorylated peptide around Y352 of human Syk (NP_003168.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - Phospho-Syk-Y352/ZAP70-Y319 pAb

Western blot analysis of extracts of Jurkat cells, using Phospho-Syk-Y352/ZAP70-Y319 antibody at 1:1000 dilution. Jurkat cells were treated by 10% FBS for 30 minutes after serum-starvation overnight. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA.