

# Phospho-Syk-Y352/ZAP70-Y319 pAb



Catalog Number: 166462

## Product name

Rabbit IgG

Phospho-Syk-Y352/ZAP70-Y319 pAb

## Applications

## Specificity

WB

Human

## Dilutions

## Antibody description


WB 1:500 - 1:2000

Polyclonal antibody to Phospho-Syk-Y352/ZAP70-Y319

## Validations

## Preparation

Antigen: A synthetic phosphorylated peptide around Y352 of human Syk (NP\_003168.2).



## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Western blot - Phospho-Syk-Y352/ZAP70-Y319 pAb

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

Western blot analysis of extracts of Jurkat cells, using Phospho-Syk-Y352/ZAP70-Y319 antibody at 1:1000 dilution. Jurkat cells were treated by 10% FBS for 30 minutes after serum-starvation overnight. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA.

## Ig Type