

Anti-ANKRD20A1 antibody



Catalog Number: 176462

Product name

Anti-ANKRD20A1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ANKRD20A1

Preparation

This antigen of this antibody was synthesized peptide derived from the internal region of human a20.

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western blot analysis of 3T3 using A20 antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000