

Anti-PGP9.5 (7B8G5E7) Mouse antibody



Catalog Number: 168462

Product name

Anti-PGP9.5 (7B8G5E7) Mouse antibody

Dilutions

WB: 1/1000

Specificity

;ICC: 1/300

Human

Antibody description

PGP9.5 (7B8G5E7) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human UCHL1 / PGP9.5 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide 50% glycerol

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG2b


Applications

WB, ICC/IF

Validations



Immunocytochemistry stain of COS7 using UCHL1 / PGP9.5 mouse mAb (1:300).



Western blot detection of UCHL1 / PGP9.5 in U251, Raji and IMR32 cell lysates and using UCHL1 / PGP9.5 mouse mAb (1:1000 diluted). Predicted band size: 25KDa. Observed band size: 25KDa.