

GANP antibody



Catalog Number: 110862

Product name

GANP antibody

IP: 1:200-1:2000

IHC: 1:10-1:100

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:20-1:200

Antibody description

GANP Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human colon cancer tissue. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 75-80 kDa

Validations




Preparation

This antibody was obtained by immunization of GANP recombinant protein (Accession Number: BC013285). Purification method: Antigen affinity purified.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110862(GANP antibody) at dilution of 1:300

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:

WB: 1:200-1:1000


Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:110862(GANP antibody) at dilution of 1:10 (under 10x lens)



GANP antibody

Catalog Number: 110862

Immunofluorescent analysis of HepG2 cells, using MCM3AP antibody Catalog No:110862 at 1:100 dilution and FITC-labeled donkey anti-rabbit IgG(green). Blue pseudocolor = DAPI (fluorescent



IP Result of anti-GANP (IP:Catalog No:110862, 4ug; Detection:Catalog No:110862 1:500) with mouse brain tissue lysate 4000ug.