

Integrin beta-1 antibody



Catalog Number: 111862

Product name

Integrin beta-1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Integrin beta-1 Rabbit Polyclonal antibody.
Positive WB detected in HeLa cells, HT-1080 cells.
Observed molecular weight by Western-blot: 130 kDa

Preparation

This antibody was obtained by immunization of Integrin beta-1 recombinant protein (Accession Number: NM_002211). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111862(ITGB1 antibody) at dilution of 1:300