

Anti-non-muscle Myosin IIB/MYH10 Rabbit antibody



Catalog Number: 167862

Product name

Anti-non-muscle Myosin IIB/MYH10 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

non-muscle Myosin IIB/MYH10 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human non-muscle Myosin IIB/MYH10

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/200

ICC/IF: 1/20

FC: 1/50

IP: 1/20

Validations



Western blot detection of non-muscle Myosin IIB/MYH10 in A549,HL-60,U251,U87-MG cell lysates using non-muscle Myosin IIB/MYH10 Rabbit mAb(1:1000 diluted).Predicted band size:229kDa.Observed band size:229kDa.