

Anti-Growth hormone receptor Rabbit antibody



Catalog Number: 168862

Product name

Rabbit IgG

Anti-Growth hormone receptor Rabbit antibody

Applications

Specificity

WB, FC, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000-1/5000

Growth hormone receptor Rabbit monoclonal antibody

FC: 1/50

IP: 1/50

Preparation

Antigen: A synthetic peptide of human Growth hormone receptor

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Western blot detection of Growth hormone receptor in Rat Brain lysates using Growth hormone receptor Rabbit mAb(1:500 diluted).Predicted band size:72kDa.Observed band size:120kDa.

