

# Anti-OR1D5 antibody



Catalog Number: 176162

## Product name

Anti-OR1D5 antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

WB, IF

## Antibody description

Rabbit polyclonal antibody to OR1D5

## Dilutions

WB: 1:500-1:2000

## Preparation

IF: 1:100-1:500

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or1d4/5. aa range:200-249

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Validations



Fig1:: Western Blot analysis of various cells using Olfactory receptor 1D4/1D5 Polyclonal Antibody diluted at 1: 2000