

TOAG1 antibody

Catalog Number: 116162

Product name

TOAG1 antibody

Specificity

Human; other species not tested.

Antibody description

TOAG1 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue, human oesophagus cancer tissue. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 35-41 kDa

Preparation

This antibody was obtained by immunization of TOAG1 recombinant protein (Accession Number: NM_173826). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116162(C3orf23 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human oesophagus cancer tissue slide using Catalog No:116162(C3orf23 Antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemistry of paraffin-embedded human oesophagus cancer tissue slide using Catalog No:116162(C3orf23 Antibody) at dilution of 1:50 (under 10x lens)


TOAG1 antibody

Catalog Number: 116162

of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:116162(C3orf23 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:116162(C3orf23 Antibody) at dilution of 1:50 (under 40x lens)