

# Anti-Bag3 Rabbit antibody

Catalog Number: 167162

## Product name

Anti-Bag3 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Bag3 Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human Bag3

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, FC

## Dilutions

WB: 1/1000

FC: 1/20

## Validations



Western blot detection of Bag3 in K562,Rat Brain cell lysates using Bag3 Rabbit pAb(1:1000 diluted).Predicted band size:62kDa.Observed band size:80kDa.