

ARP3 antibody

Catalog Number: 108162

Product name

ARP3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ARP3 Rabbit Polyclonal antibody. Positive IP detected in mouse spleen tissue. Positive WB detected in mouse spleen tissue, human liver tissue. Observed molecular weight by Western-blot: 47kd

Preparation

This antibody was obtained by immunization of ARP3 recombinant protein (Accession Number: NM_005721). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:108162(ARP3 antibody) at dilution of 1:1000



IP Result of anti-ARP3 (IP:Catalog No:108162, 4ug; Detection:Catalog No:108162 1:1000) with mouse spleen tissue lysate 4000ug.