

MCSF antibody



Catalog Number: 112562

Product name

MCSF antibody

Specificity

human; other species not tested.

Antibody description

MCSF Rabbit Polyclonal antibody. Positive WB detected in U-937 cells. Observed molecular weight by Western-blot: 70kd

Preparation

This antibody was obtained by immunization of MCSF recombinant protein (Accession Number: NM_172212). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



U-937 cells were subjected to SDS PAGE followed by western blot with Catalog No:112562(MCSF antibody) at dilution of 1:200