

Anti-Phospho-Rb2 p130 (Thr986) Rabbit antibody



Catalog Number: 168562

Product name

Anti-Phospho-Rb2 p130 (Thr986) Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse

Applications

WB, FC, IP

Antibody description

Phospho-Rb2 p130 (Thr986) Rabbit monoclonal antibody

Dilutions

1:1000

1:50

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Thr986 of human p130

1:20

Validations

Formulation


Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal



Western blot detection of Phospho-Rb2 p130 (Thr986) in K562 cell lysates using Phospho-Rb2 p130 (Thr986) Rabbit mAb(1:1000 diluted). Predicted band size:128kDa. Observed band size:130kDa.