

SYTL4 antibody

Catalog Number: 115962

Product name

SYTL4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SYTL4 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of SYTL4 recombinant protein (Accession Number: NM_001370167). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115962(SYTL4 Antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:115962(SYTL4 antibody) at dilution of 1:50 (under 10x lens)