

# ZNF141 Polyclonal Antibody



Catalog Number: 162262

---

## Product name

ZNF141 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to ZNF141

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 405-474 of human ZNF141 (NP\_003432.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


IHC, IF

## Dilutions

IHC 1:50 - 1:200

IF 1:50 - 1:200

## Validations



Immunohistochemistry - ZNF141 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using ZNF141 antibody at dilution of 1:100 (40x lens).