

HEXIM1 antibody



Catalog Number: 107262

Product name

HEXIM1 antibody

Specificity

Human; other species not tested.

Antibody description

HEXIM1 Mouse Monoclonal antibody. Positive WB detected in human testis tissue, MCF-7 cells. Observed molecular weight by Western-blot: 54 kDa

Preparation

This antibody was obtained by immunization of HEXIM1 recombinant protein (Accession Number: NM_006460). Purification method: Thiophilic affinity chromatograph.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgM

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107262(HEXIM1 Antibody) at dilution of 1:1000