

# Anti-KCNIP3 antibody



Catalog Number: 142662

---

## Product name

Anti-KCNIP3 antibody

## Immunogen

[Human CSEN / calsenilin / KCNIP3 \(His Tag\) recombinant protein](#)

## Specificity

Human CSEN / calsenilin / KCNIP3

## Antibody description

Rabbit Monoclonal to Human KCNIP3

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human CSEN / calsenilin / KCNIP3 .

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 1-5 µg/mL

## Validations



Human CSEN / calsenilin / KCNIP3  
Immunohistochemistry(IHC) 14555