

Anti-PRIMPOL antibody



Catalog Number: 176662

Product name

Anti-PRIMPOL antibody

Specificity

Human

Antibody description

Mouse monoclonal antibody to PRIMPOL

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ICC, IHC-P

Dilutions

ICC: 1:100-1:200

IHC-P: 1:200

Validations




Fig1: ICC staining CCDC111 in HepG2 cells(red). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Fig2: Immunohistochemical analysis of paraffin-embedded human colonic tissue using anti-CCDC111 mouse monoclonal antibody.