

# UBR1 antibody

Catalog Number: 116662

## Product name

UBR1 antibody

## Specificity

Human; other species not tested.

## Antibody description

UBR1 Rabbit Polyclonal antibody. Positive IP detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 200 kDa

## Preparation

This antibody was obtained by immunization of UBR1 recombinant protein (Accession Number: NM\_174916). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:200-1:1000

## Validations



IP Result of anti-UBR1 (IP:Catalog No:116662, 4ug; Detection:Catalog No:116662 1:300) with mouse skeletal muscle tissue lysate 4000ug.