

C14orf2 antibody

Catalog Number: 108662

Product name

C14orf2 antibody

Specificity

Human; other species not tested.

Antibody description

C14orf2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, human brain tissue. Positive IHC detected in human heart tissue. Observed molecular weight by Western-blot: 7kd

Preparation

This antibody was obtained by immunization of C14orf2 recombinant protein (Accession Number: NM_004894). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108662(C14orf2 antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:108662(C14orf2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:108662(C14orf2 Antibody) at dilution of 1:50 (under 40x lens)