

Mouse KLK11/Kallikrein 11 (His Tag) recombinant protein



Catalog Number: 500362

General Information

Gene Name Synonym

Hippostasin; Serine protease 20

Protein Construction

A DNA sequence encoding the mouse KLK11 isoform 2 (Q9QYN3) (Met 1-Asn 276), was fused with a C-terminal polyhistidine tag.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu 45

Molecular Mass

The secreted recombinant mouse KLK11 (pro

form) comprises 243 amino acids and has a calculated molecular mass of 27.3 kDa. The apparent molecular mass of rmKLK11 is approximately 38-42 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse KLK11 / Kallikrein-11 Protein (His Tag)
SDS-PAGE