

# HIST1H4F antibody



Catalog Number: 111362

## Product name

HIST1H4F antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HIST1H4F Rabbit Polyclonal antibody. Positive WB detected in human spleen tissue. Observed molecular weight by Western-blot: 11 kDa

## Preparation

This antibody was obtained by immunization of HIST1H4F recombinant protein (Accession Number: NM\_003540). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:111362(HIST1H4F antibody) at dilution of 1:300