

# KIFAP3 antibody



Catalog Number: 107362

## Product name

KIFAP3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

KIFAP3 Mouse Monoclonal antibody. Positive WB detected in human testis tissue. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 91-100kd

## Preparation

This antibody was obtained by immunization of KIFAP3 recombinant protein (Accession Number: NM\_014970). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107362(KIFAP3 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:107362(KIFAP3 Antibody) at dilution of 1:50