

Anti-Cytokeratin 19 (4E8) Mouse antibody



Catalog Number: 168362

Product name

Anti-Cytokeratin 19 (4E8) Mouse antibody

Specificity

Human

Antibody description

Cytokeratin 19 (4E8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human KRT19 expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IHC-P, ICC/IF

Dilutions


WB: 1/500 - 1/2000

IHC: 1/200 - 1/1000

ICC: 1/200 - 1/1000

ELISA: 1/10000

Validations



Western blot detection of Cytokeratin 19 in CaCo2, HepG2 and T-47D cell lysates using Cytokeratin 19 mouse mAb (dilution 1:500). Predicted band size: 41KDa. Observed band size: 41KDa.