

# MT3 antibody

Catalog Number: 112762

---

## Product name

MT3 antibody

## Immunogen

[Human MT3 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MT3 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue, human kidney tissue, human testis tissue.

## Preparation

This antibody was obtained by immunization of MT3 recombinant protein (Accession Number: NM\_005954). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, IHC


## Dilutions

Recommended Dilution:  
IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human brain using Catalog No:112762(MT3 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:112762(MT3 antibody) at dilution of 1:50 (under 40x lens)