

# PDAP1 antibody

Catalog Number: 113762

## Product name

PDAP1 antibody

IF: 1:10-1:100


## Specificity

human, mouse, rat; other species not tested.

## Antibody description

PDAP1 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells. Positive IHC detected in human breast cancer. Positive IF detected in Hela. Observed molecular weight by Western-blot: 28kd

## Validations




## Preparation

This antibody was obtained by immunization of PDAP1 recombinant protein (Accession Number: NM\_014891). Purification method: Antigen affinity purified.

MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:113762(PDAP1 antibody) at dilution of 1:800

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT


## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:113762(PDAP1 antibody) at dilution of 1:50 (under 10x lens)

## Ig Type

Rabbit IgG



## Applications

ELISA, WB, IHC, IF

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200


Immunohistochemical of paraffin-embedded human breast cancer using Catalog

# PDAP1 antibody

Catalog Number: 113762

---

No:113762(PDAP1 antibody) at dilution of 1:50  
(under 40x lens)



Immunofluorescent analysis of HeLa cells, using PDAP1 antibody Catalog No:113762 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).