

Anti-Visfatin Rabbit antibody



Catalog Number: 167762

Product name

Anti-Visfatin Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Visfatin Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Visfatin

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/200

ICC/IF: 1/20

FC: 1/50

IP: 1/20

Validations



Western blot detection of Visfatin in HeLa, A549, HL-60, U2OS, C6 cell lysates using Visfatin Rabbit mAb(1:1000 diluted). Predicted band size: 56kDa. Observed band size: 56kDa.