

Anti-Importin 9/RANBP9 Rabbit antibody



Catalog Number: 168762

Product name

Rabbit IgG

Anti-Importin 9/RANBP9 Rabbit antibody

Applications

Specificity

WB, ICC/IF, FC, IP

Human, Mouse

Dilutions

Antibody description

WB: 1/1000

Importin 9/RANBP9 Rabbit monoclonal antibody

ICC/IF: 1/50

Preparation

FC: 1/50

Antigen: A synthetic peptide of human Importin 9

IP: 1/20


Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Validations

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.



Clonality

Monoclonal

Western blot detection of Importin 9/RANBP9 in Jurkat, HeLa cell lysates using Importin 9/RANBP9 Rabbit mAb(1:1000 diluted). Predicted band size: 115kDa. Observed band size: 115kDa.

Ig Type