

Human DPP3 (His tag) recombinant protein



Catalog Number: 519762

General Information

Gene Name Synonym

Dipeptidyl peptidase 3, Dipeptidyl aminopeptidase III, Dipeptidyl arylamidase III, Dipeptidyl peptidase III, Enkephalinase B, DPP3

Protein Construction

A DNA sequence encoding the Human DPP3 Met1-Ala737 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human DPP3 (His tag) recombinant protein