

# Mouse REG3D (His Tag) recombinant protein



Catalog Number: 500072

## General Information

### Gene Name Synonym

Reg III delta; Reg3d protein; Regenerating islet-derived 3 delta, isoform CRA\_a

### Protein Construction

A DNA sequence encoding the mouse REG3D isoform 1 (NP\_038921.2) extracellular domain (Met 1-Gly 175) was fused with a polyhistidine tag at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

> 97 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Glu 27

### Molecular Mass

The secreted recombinant mouse REG3D comprises 160 amino acids and has a predicted molecular mass of 18.5 kDa. rmREG3D migrates as an approximately 17 kDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Mouse REG3D Protein (His Tag) SDS-PAGE