

# GPR107 antibody



Catalog Number: 111072

## Product name

GPR107 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

GPR107 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, mouse cerebellum tissue. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 60-67 kDa

## Preparation

This antibody was obtained by immunization of GPR107 recombinant protein (Accession Number: NM\_020960). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111072(GPR107 Antibody) at dilution of 1:600



Immunofluorescent analysis of HEK-293 cells using Catalog No:111072(GPR107 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG