

# NDUFB2 antibody



Catalog Number: 113072

## Product name

NDUFB2 antibody

Recommended Dilution:

WB: 1:500-1:5000

## Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

## Antibody description

NDUFB2 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue, human heart tissue, human kidney tissue, human liver tissue, human ovary tissue, human placenta tissue, human skin tissue, human spleen tissue, human testis tissue. Positive WB detected in human heart tissue. Observed molecular weight by Western-blot: 12kd

## Validations




## Preparation

This antibody was obtained by immunization of NDUFB2 recombinant protein (Accession Number: BC063026). Purification method: Antigen affinity purified.

human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113072(NDUFB2 antibody) at dilution of

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



1:1000

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human brain using Catalog No:113072(NDUFB2 antibody) at dilution of 1:100 (under 10x lens)


## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions



# NDUFB2 antibody



Catalog Number: 113072

---

Immunohistochemical of paraffin-embedded

human brain using Catalog No:113072(NDUFB2 antibody) at dilution of 1:100 (under 40x lens)