

# COL10A1 Polyclonal Antibody



Catalog Number: 165072

## Product name

COL10A1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to COL10A1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 521-680 of human COL10A1 (NP\_000484.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC

## Dilutions


IHC 1:100 - 1:200

## Validations



Immunohistochemistry - COL10A1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat lung using COL10A1 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - COL10A1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver using COL10A1 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - COL10A1 Polyclonal Antibody


# COL10A1 Polyclonal Antibody



Catalog Number: 165072

---

Immunohistochemistry of paraffin-embedded human stomach using COL10A1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - COL10A1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using COL10A1 antibody at dilution of 1:100 (40x lens).