

SDR42E1 antibody

Catalog Number: 115072

Product name

SDR42E1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SDR42E1 Rabbit Polyclonal antibody. Positive IHC detected in human bladder tissue.

Preparation

This antibody was obtained by immunization of SDR42E1 recombinant protein (Accession Number: BC019670). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human bladder tissue slide using Catalog No:115072(SDR42E1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human bladder tissue slide using Catalog No:115072(SDR42E1 Antibody) at dilution of 1:50 (under 40x lens)