

ATP5C1 antibody



Catalog Number: 107072

Product name

ATP5C1 antibody

Immunogen

[Human ATP5C1 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ATP5C1 Mouse Monoclonal antibody. Positive WB detected in human heart tissue, fetal human brain tissue. Observed molecular weight by Western-blot: 33 kDa

Preparation

This antibody was obtained by immunization of ATP5C1 recombinant protein (Accession Number: XM_047425254). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:107072(ATP5C1 Antibody) at dilution of 1:1000