

Anti-N Cadherin Rabbit antibody



Catalog Number: 169072

Product name

Anti-N Cadherin Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

N Cadherin Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human N Cadherin

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions


WB: 1/500-1/2000

IHC: 1/50-1/200

Validations



Western blot analysis of N-Cadherin expression in



HeLa cell lysate.

Immunohistochemical analysis of paraffin-embedded human colon, using N Cadherin Antibody.