

FAM133A antibody



Catalog Number: 110472

Product name

FAM133A antibody

Immunogen

[Human FAM133A Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

FAM133A Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue.

Preparation

This antibody was obtained by immunization of FAM133A recombinant protein (Accession Number: XM_011530926). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:110472(FAM133A Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:110472(FAM133A Antibody) at dilution of 1:50 (under 40x lens)