

Anti-PLGRKT antibody



Catalog Number: 176472

Product name

Polyclonal

Anti-PLGRKT antibody

Ig Type

Specificity

Rabbit IgG

Human, Rat

Applications

Antibody description

WB

Rabbit polyclonal antibody to PLGRKT

Dilutions

Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human plgrkt. aa range:1-50

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Validations




Fig1.; Western Blot analysis of HEPG2 cells using Plasminogen receptor Polyclonal Antibody diluted at 1:500. Secondary antibody (catalog#: HA1001) was diluted at 1:20000