

Anti-PLGRKT antibody

Catalog Number: 176472

Product name

Anti-PLGRKT antibody

Specificity

Human, Rat

Antibody description

Rabbit polyclonal antibody to PLGRKT

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human plgrkt. aa range:1-50

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of HEPG2 cells using Plasminogen receptor Polyclonal Antibody diluted at 1:500. Secondary antibody (catalog#: HA1001) was diluted at 1:20000