

# UBAC2 antibody

Catalog Number: 116472

## Product name

UBAC2 antibody

## Immunogen

[Human UBAC2 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

UBAC2 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, HEK-293 cells, mouse kidney tissue. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 35 kDa

## Preparation

This antibody was obtained by immunization of UBAC2 recombinant protein (Accession Number: XM\_047430286). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116472(UBAC2 Antibody) at dilution of 1:600



Immunofluorescent analysis of HEK-293 cells using Catalog No:116472(UBAC2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)